


eliq step




Load alert orthopedic device:

- Conservative treatment
- Postoperative treatment



Why use Cliqstep?

- Cliqstep is a overload, sound alerting device that gives you a personal rehabilitation training with different load resistance.
- Several orthopedic condition (prostheses, fractures, muscle and rupture of tendons) of the lower extremities requires gradually load of the foot (11 lbs., 22 lbs., 33 lbs.) over number of weeks to increase the leg strength. By using Cliqstep there will be less risk of overload the injury during rehabilitation period.
- Cliqstep is made by recyclable material.



How to attach Cliqstep?

- Cliqstep can be mounted on plaster foot by wrapping the bandage as fixation.
- To mount Cliqstep on shoe splint remove the flanges and attach with 4 screws.



30 kg



20 kg



12 kg

Product description:

Exchangeable load sensors.

Red: 30 kg.

Yellow: 20 kg.

Blue: 12 kg.

Dark grey holder

Flat spring: generate “click” sound when exposed to pressure

eliq step



Customer feedback from user tests:

Orthopedists: "Easy to attach"

Physiotherapist: "Works perfectly"

Physiotherapist: "Smart solution to old problem"

Patient: "No problem to adapt to this rehabilitation procedure"

Druiden AB

Jan-Olof Svärd

info@druiden.se